English	Instructions for Use: Dental Instruments	2
Español	Instrucciones de uso: Instrumentos Dentales	6
русский	Инструкция по применению: Зубные протезы	10
Deutsch	Gebrauchsanweisung: Zahnmedizinische Instrumente	14
Français	Instructions d'utilisation: Instruments Dentaires	18
Italiano	Istruzioni per l'uso: Strumenti Dentali	22
Português	Instruções de utilização: Instrumentos Dentários	26
Türkçe	Kullanım talimatları: Dental Aletler	30
简体中文	使用说明：牙科义齿	34
日本語	使用説明書：歯科プロテーゼ	37
한국어	적응증: 치과 보철	41
عربي	إرشادات الاستخدام: أدوات طب الأسنان	45
Polski	Instrukcja obsługi: Narzędzia Dentystyczne	49

BioHorizons
2300 Riverchase Center
Birmingham, AL 35244 USA
TOLL FREE 888.246.8338
TEL 205.967.7880 FAX 205.870.0304
www.biohorizons.com

Altatec GmbH
Maybachstrasse 5
71299 Wimsheim
Germany
TEL +49.7044.94450

2797



*The CE mark is valid
only if it is also printed
on the product label.*

L02069 Rev B AUG 2021

English

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic Instructions for Use
	Manufacturer
	Intra-Lock products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive (93/43/EEC) as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745. The CE mark is valid only if it is also printed on the product label.
	Reference/article number
	Lot/batch number
	Do not re-use
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
	Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.
	EU Authorized Representative
	Do not use if package is damaged
	Medical Device
	Non-sterile
	Do not re-sterilize
	Home position

This document supersedes all prior revisions. Original language is English.

This document applies to Intra-Lock Dental Instruments and Accessories. Every product packaging label contains a description of the enclosed device.

DESCRIPTION

Intra-Lock instruments are used for dental implant procedures such as site development, implant placement, and implant restorations within the specific indications of each implant system. The label on each instrument contains important product information including whether the instrument is supplied sterile or non-sterile and if the instrument is for single-use only. Instruments and kits supplied non-sterile must be cleaned and sterilized prior to first use and each use thereafter.

INDICATIONS FOR USE

Intra-Lock instruments and surgical kits are indicated for use in site development, placement, and restoration of Intra-Lock implants and associated components.

CONTRAINDICATIONS

Intra-Lock instruments should not be used with patients that have allergies to the specific materials used, including stainless steel and titanium alloy.

DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Intra-Lock strongly recommends completion of dental implant courses and strict adherence to the instructions pertaining to Intra-Lock products.

Drills should be replaced when wear is noticed, such as a decrease in cutting efficiency or when signs of discoloration appear. Cutting instruments should be replaced after approximately 12 to 20 osteotomy cycles, depending on bone density. Intra-Lock recommends use of a drill usage chart to track drill use and to ensure drills are replaced as directed. Guided Drills and Drill Guides: Care must be taken to ensure the Guided Drills follow the angle of the guide handle inserts or sleeves to prevent drill binding. Refer to the drill report for proper drill selection and usage.

WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any Intra-Lock Instructions for Use (IFU). Refer to OSHA standard 29CFR1910.1030 prior to cleaning and sterilization. Additional technical information is available upon request from Intra-Lock, or may be viewed and/or downloaded at www.intra-lock.com. Contact Intra-Lock Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

Prosthetic Instruments - Only use a torque wrench for final installation of prosthetic components requiring a specified installation torque. Finger-tightening of prosthetic components can result in insufficient torque and eventual loosening of the component. Over-tightening prosthetic components could break the component or spin the mating implant.

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with instruments, prosthetics, and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) abutment screw and/or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

HANDLING AND STERILIZATION

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. If the product is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined. If the product is supplied non-sterile or when reprocessing: Remove and discard any shipping material before initial sterilization. Assembled instruments (e.g. Overdenture Cap Tool) should be disassembled before each cleaning and sterilization cycle to avoid debris encapsulation, material discoloration and/or inappropriate drying of components. Disassemble the surgical kits before each cleaning. The device must be cleaned and sterilized. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for ten (10) minutes minimum.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

(Guided Surgery Drill Guide Handle: blow pressurized air through one end to purge. For both instruments: use a medical instrument lubricant such as STERIS' Hinge Free® after each use.) Refer to the labeling of the cleaning/lubricating agent used for instructions for use. Instruments that are unable to be cleaned, discolored, do not properly interface with mating components, and/or do not articulate as designed should be disposed of.

If applicable, return the instruments to the appropriate locations in the surgical tray. The Intra-Lock Prosthetics Tray should be cleaned and sterilized empty. For sterilization, place the cleaned device in an approved sterilization bag or wrap and run through one of the following qualified sterilization cycles:

Sterilization Cycles			
Reference:	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	UK HTM 01-01 Part C:2016
Type:	Gravity Steam	Prevacuum Steam	Prevacuum Steam
Exposure Time and Temperature:	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
Minimum Dry Time:	30 minutes	20 minutes	20 minutes

The following qualified sterilization cycle may also be used for individual instruments, the Intra-Lock Surgical Kit, and the Intra-Lock Surgical Tray:

Sterilization Cycles	
Reference:	ANSI/AAMI TIR12:2010
Type:	Gravity Steam
Exposure Time and Temperature:	30 minutes at 121°C (250°F)
Minimum Dry Time:	30 minutes

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

Attention! Improper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

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