<b>English</b>	Instructions for Use: Dental Implants	2
<b>Español</b>	Instrucciones de uso: implantes dentales	5
<b>Русский</b>	Инструкция по эксплуатации: зубные имплантаты	8
<b>Deutsch</b>	Gebrauchsanweisung: Zahnimplantate	12
<b>Français</b>	Instructions d'utilisation : Implants dentaires	16
<b>Italiano</b>	Istruzioni per l'uso: impianti dentali	20
<b>Português</b>	Instruções de utilização: Implantes Dentários	23
<b>Türkçe</b>	Kullanım Talimatları: Dental İmplantler	26
<b>简体中文</b>	使用说明：牙植入体	29
<b>日本語</b>	使用説明書・歯科インプラント	32
<b>한국어</b>	사용 지침: 치과 수술용 임플란트	35
<b>ع العربية</b>	إرشادات الاستخدام: غرسات الأسنان	38
<b>Polski</b>	Instukcja użycia: implanty dentystyczne	41

**ARCHIVE**

BioHorizons  
2300 Riverchase Center  
Birmingham, AL 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
[www.biohorizons.com](http://www.biohorizons.com)

Altatec GmbH  
Maybachstrasse 5  
71299 Wimsheim  
Germany  
TEL +49.7044.94450

2797

*The CE mark is valid  
only if it is also printed  
on the product label.*

L01079 Rev B OCT 2022

## English

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic Instructions for Use
	Manufacturer
	Intra-Lock products may carry the CE mark and fulfill the requirements of the Medical Devices Directive 93/42/EEC. The CE mark is valid only if it is also printed on the product label.
	Reference/article number
	Lot/batch number
	Do not re-use
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution, and use by, or on the order of, a dentist or physician.
	EU Authorised Representative
	Do not use if package is damaged
	Medical Device
<b>Non-sterile</b>	Non-sterile
	Do not re-sterilize

This document supersedes all prior revisions. Original language is English.

This document applies to Intra-Lock Dental Implants.

## **DESCRIPTION**

Intra-Lock implants are manufactured from biocompatible titanium alloy and are single use only.

## **INDICATIONS FOR USE**

Intra-Lock implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion or (2) when splinted together for multiple tooth replacement or when stabilized with an overdenture supported by multiple implants.

Intra-Lock 3.3mm diameter implants may be used as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion, (2) when splinted together as an artificial root structure for multiple tooth replacement of mandibular incisors, or (3) for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be placed in immediate function when good primary stability has been achieved and with appropriate occlusal loading.

## **CONTRAINdications**

Intra-Lock implants should not be used in patients who have contraindicating systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity.

Intra-Lock 3.3mm diameter implants are contraindicated for single tooth replacement of maxillary central incisors.

## **DIRECTIONS FOR USE**

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. Intra-Lock strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany Intra-Lock implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinse on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55:25-30, 1997, Suppl5).

**IMPORTANT!** (1) The implant procedure manual should be reviewed prior to attempting placement. (2) In dense (D1-D2) bone, further preparation of the site may be required to avoid applying excessive torque to the implant connection. (3) To further prepare the site, the crestal bone drills may be used. (4) The clinician should verify that the proper instruments are available and at hand prior to attempting placement.

## **WARNINGS AND PRECAUTIONS**

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone width and depth are available at the desired implant site and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth.

Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments. For short implants, clinicians should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting the short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration avoid immediate loading. Intra-Lock dental implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of Intra-Lock dental implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but not limited to: casting beyond a 30-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. Small diameter implants with either angled or straight abutments are intended for the anterior region of the mouth and not intended for the posterior region of the mouth due to possible failure of the implant.

Devices should not be used in patients with a known allergy or sensitivity to the device material.

#### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histological responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

#### **HANDLING AND STERILIZATION**

This product is supplied sterile and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. Implants are single use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. Intra-Lock assumes no responsibility for attempted re-use or re-sterilization.

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